



## Face Mask Nose Clip

 Vratik

[VIEW IN BROWSER](#)

updated 24. 4. 2020 | published 24. 4. 2020

### Summary

Nose clip for face mask. Better isolation around the nose and prevent of fogging glasses.

[Healthcare](#) > [Accessibility](#)

covidmask

covid19protection

covid19

Nose clip for face mask. Better isolation around the nose and prevent of fogging glasses.

If you dont find correct size, leave a comment.

Version 1.1 - longest "teeths" with small bead for better fit

Version 1.3 - wider front tooth; radiuses added

Version 1.3\_small - for smaller nose :-)

Version 1.3\_wide - for wider nose

FLAT Version added:

Printed from PLA can be shaped using hot water:

1. Print from PLA! (100%infill)
2. Prepare some container (plate, bowl ...), paper towel and hot water.
3. Pour hot water over printed part and wait few seconds (10 or 15).
4. Carefully pull out part from water (use tweezers) and dry it quickly (to avoid burning by hot water.)
5. Place it over your nose (down worry, dry part wont burn you) and shape it.
6. Wait until it cool down (1 minute) or use cold wet napkin.
7. Enjoy your personal nose clip :-)

## Model Files (.stl, .3mf, .obj, .amf)

[↓ DOWNLOAD ALL FILES](#)

**facemasknoseclip13.stl**

updated 24. 4. 2020

853.5 KB



**facemasknoseclip13\_small.stl**

updated 24. 4. 2020

860.7 KB



**facemasknoseclip13\_wide.stl**

updated 24. 4. 2020

689.7 KB



**noseclipflat2mm.stl**

updated 24. 4. 2020

426.6 KB



**noseclipflat38mm.stl**

updated 24. 4. 2020

399.0 KB



**noseclipflat3mm.stl**

updated 24. 4. 2020

399.5 KB





The Author has not uploaded any print files.  
Try to search in [User print files](#) section or generate and upload your own.

## License

This work is licensed under a  
**Creative Commons (4.0 International License)**



**Attribution-NonCommercial**

- 
- ✗ | Sharing without ATTRIBUTION
  - ✓ | Remix Culture allowed
  - ✗ | Commercial Use
  - ✗ | Free Cultural Works
  - ✗ | Meets Open Definition